

# Quantitative Analysis of 25-OH-D<sub>2</sub> and 25-OH-D<sub>3</sub> in Serum by LC-MS/MS

Marta Kozak, Thermo Fisher Scientific, San Jose, CA

## Key Words

- TSQ Quantum Access™
- Accela™ LC System
- Quantitation
- Vitamin D

## Introduction

Vitamin D is a group of fat-soluble prohormones that helps the body absorb calcium and maintain proper levels of calcium and phosphorus in the blood. A vitamin D deficiency results in impaired bone mineralization and leads to bone softening diseases, such as rickets in children and osteomalacia and osteoporosis in adults, as well as autoimmune and nervous system diseases. A vitamin D overdose can raise blood levels of calcium and cause gastrointestinal symptoms and kidney disease. High blood levels of calcium also can cause heart rhythm abnormalities and calcinosis, the deposition of calcium and phosphate in the body's soft tissues.

Vitamin D exists naturally in two major forms: vitamin D<sub>2</sub> and vitamin D<sub>3</sub>. (See Figure 1.) Vitamin D<sub>2</sub> is derived from plant and fungal sources. Vitamin D<sub>3</sub> is derived from animal sources and is produced in skin exposed to ultraviolet (UVB) radiation from the sun. In the body, vitamin D is metabolized to 25-hydroxy-vitamin D (25-OH-D) in the liver and subsequently to biologically active 1,25-dihydroxy-vitamin D (1,25(OH)<sub>2</sub>D) in the kidneys.

The level of vitamin D in serum is measured by analyzing its major metabolites. Among the vitamin D metabolites, 25-OH-D has the highest concentration and longest half-life in blood. It is considered to be the best indicator of vitamin D concentration in the blood. Serum 25-OH-D concentrations in the range 32–80 ng/mL (80–200 nMol/L) are optimum. Concentrations below 32 ng/mL indicate a vitamin D deficiency, and concentrations above 100 ng/mL are regarded as toxic.

## Goal

To develop a quantitative, efficient LC-MS method for analysis of 25-hydroxy-vitamin D<sub>3</sub> and 25-hydroxy-vitamin D<sub>2</sub> in serum.




Figure 1: 25-hydroxy-vitamin-D<sub>3</sub> and 25-hydroxy-vitamin-D<sub>2</sub> chemical structures.

## Experimental Conditions

Figure 2 illustrates the workflow for the quantitative analysis of 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> in serum.

### Chemicals and Reagents

The reagents 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> were purchased from Sigma-Aldrich (Sigma-Aldrich, St. Louis, MO) and were stored at –20°C. The internal standard [<sup>2</sup>H<sub>6</sub>]-25-OH-D<sub>3</sub> was purchased from Medical Isotopes (Medical Isotopes, Inc., Pelham, NH) and was stored at –70°C. Bovine serum albumin and PBS buffer were purchased from Sigma-Aldrich and were stored in a refrigerator. Bovine serum was used because human serum with low 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> content was not commercially available.

### Sample Preparation

To prepare each sample for analysis, 100 µL of internal standard solution (100 ng/mL [<sup>2</sup>H<sub>6</sub>]-25-OH-D<sub>3</sub> in ethanol) and 1 mL of acetonitrile were added to 250 µL of serum. The sample was vortexed and centrifuged. The supernatant was decanted and evaporated to dryness at 37°C. The sample was then reconstituted in 100 µL of methanol. After reconstitution, 25 µL was injected into the LC-MS.

### Calibration Curve Standards Preparation

A standard stock solution of 100 µg/mL 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> in ethanol was prepared. Standard spiking solutions of 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> in ethanol at concentrations of 10 µg/mL, 1 µg/mL, and 100 ng/mL were prepared by serial dilution. The appropriate amount of the standard spiking solution was added to 250 µL of 6% BSA in 0.01M PBS (pH 7.4) to prepare calibration standards at the following eight concentrations: 4 ng/mL, 10 ng/mL, 20 ng/mL, 40 ng/mL, 100 ng/mL, 200 ng/mL, 400 ng/mL, and 500 ng/mL. The standards were processed




Figure 2: Workflow for the quantitative analysis of 25-OH-D<sub>3</sub> and 25-OH-D<sub>2</sub> in serum.

with the sample preparation procedure described above. The standard stock solution and the standard spiking solutions were stored at  $-20^{\circ}\text{C}$ .

### HPLC

HPLC analysis was performed using the Accela High Speed LC System (Thermo Fisher Scientific, San Jose, CA). The 25  $\mu\text{L}$  samples were injected directly onto a 2.1 mm ID $\times$ 50 mm column, packed with 5  $\mu\text{m}$ , C18-aQ stationary phase. Mobile phase A was water containing 0.1% formic acid. Mobile phase B was methanol containing 0.1% formic acid. The flow rate was 400  $\mu\text{L}/\text{min}$ . The gradient is described in Table 1.

Time (min)	% A	% B
0	40	60
0.5	40	60
1.5	2	98
2.5	2	98
2.6	40	60
3.5	40	60

Table 1: Gradient profile

### MS

MS analysis was carried out on a TSQ Quantum Access triple stage quadrupole mass spectrometer with an Ion Max<sup>™</sup> ion source and an atmospheric pressure chemical ionization (APCI) probe (Thermo Scientific, San Jose, CA). The MS conditions were as follows:

Ion polarity: Positive ion mode  
Vaporizer temperature:  $375^{\circ}\text{C}$   
Capillary temperature:  $250^{\circ}\text{C}$   
Discharge current: 2.0  $\mu\text{A}$   
Sheath gas pressure ( $\text{N}_2$ ): 35 units  
Auxiliary gas pressure ( $\text{N}_2$ ): 5 units  
Scan type: H-SRM  
Scan time: 0.050 s

Table 2 lists the SRM transitions and their parameters. The parent molecular ions for the three analytes show a significant ion due to the loss of a water molecule from the molecular ion. This allows the use of the MH-18  $m/z$  peak for quantitation by SRM.

Analyte	Parent Mass	Product Mass	Collision Energy	Tube Lens
25-OH- $\text{D}_3$	383.28*	211.094	21	104
25-OH- $\text{D}_2$	395.29*	209.039	23	102
25-OH- $\text{D}_3$ (Internal Standard)	389.341*	211.157	31	103

\*parent molecular ion minus water [MH-18]

Table 2: SRM transitions

## Results and Discussion

Representative - SRM chromatograms for 25-OH- $\text{D}_3$  and 25-OH- $\text{D}_2$  at 4 ng/mL and 100 ng/mL are shown in Figures 3 and 4, respectively. Clearly identifiable and quantifiable peaks were observed.

Figures 5 and 6 show the linear fit calibration curves for 25-OH- $\text{D}_3$  and 25-OH- $\text{D}_2$ , respectively. The calibration curves have  $R^2$  values greater than 0.999, which indicate excellent linear fits over the dynamic range of 4–500 ng/mL for 25-OH- $\text{D}_3$  and 2–500 ng/mL for 25-OH- $\text{D}_2$ . The LOQ value is 4 ng/mL for 25-OH- $\text{D}_3$  and 2 ng/mL for 25-OH- $\text{D}_2$ , with LOD values  $\sim 4$  times lower for each analyte.




Figure 3: Chromatograms of 4 ng/mL standard




Figure 4: Chromatograms of 100 ng/mL standard




Figure 5: 25-OH-D<sub>3</sub> calibration curve




Figure 6: 25-OH-D<sub>2</sub> calibration curve

The method accuracy and precision were evaluated by analyzing quality control samples prepared at concentrations of 10 ng/mL, 100 ng/mL, and 400 ng/mL in 6% BSA. Intra-assay variability was determined by processing and analyzing five replicates of each QC sample. Inter-assay variability was determined by processing and analyzing five replicates of each QC sample in three different batches. Intra-assay and inter-assay variability are <8.7% and <5.8% for 25-OH-D<sub>3</sub>, and <3.2% and <2.5% for 25-OH-D<sub>2</sub>, respectively, meeting industry criteria.

## Conclusion

A fast, sensitive and reliable LC-MS/MS-SRM method has been developed for the determination of 25-OH D<sub>2</sub> and 25-OH-D<sub>3</sub> in serum. Sample analysis was performed with a runtime of less than four minutes with a quantification limit of 4 ng/mL (10 nMol/L) for 25-OH-D<sub>3</sub>, 2 ng/mL (4.85 nMol/L) for 25-OH-D<sub>2</sub> and a linearity range of 4–500 ng/mL (10–1250 nMol/L) for 25-OH D<sub>3</sub>, 2–500 ng/mL (4.85–1212.7 nMol/L) for 25-OH-D<sub>2</sub>. Data was acquired using H-SRM. The low intra-assay and inter-assay variability of the results demonstrates the reliability of the method.

## Legal Notices

©2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**View additional Thermo Scientific LC/MS application notes at: [www.thermo.com/appnotes](http://www.thermo.com/appnotes)**

*In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.*

**Africa**  
+43 1 333 5034 127  
**Australia**  
+61 2 8844 9500  
**Austria**  
+43 1 333 50340  
**Belgium**  
+32 2 482 30 30  
**Canada**  
+1 800 530 8447  
**China**  
+86 10 5850 3588  
**Denmark**  
+45 70 23 62 60  
**Europe-Other**  
+43 1 333 5034 127  
**France**  
+33 1 60 92 48 00  
**Germany**  
+49 6103 408 1014  
**India**  
+91 22 6742 9434  
**Italy**  
+39 02 950 591  
**Japan**  
+81 45 453 9100  
**Latin America**  
+1 608 276 5659  
**Middle East**  
+43 1 333 5034 127  
**Netherlands**  
+31 76 587 98 88  
**South Africa**  
+27 11 570 1840  
**Spain**  
+34 914 845 965  
**Sweden/Norway/Finland**  
+46 8 556 468 00  
**Switzerland**  
+41 61 48784 00  
**UK**  
+44 1442 233555  
**USA**  
+1 800 532 4752

[www.thermo.com](http://www.thermo.com)



Thermo Fisher Scientific,  
San Jose, CA USA is ISO Certified.

AN62452\_E 10/07S